

Work happily




V1.0


naggura'swop

naggura.com

Our products


 **Podiatry**

S503WHITE S503GREY




The Swop 503 is a 3-engine podiatry bed of optimum stability and ultra-compact frame.

Wood, Metal or Pmma




Collections


Nordic
METAL



Yellow



Blue




Pink

Foll
METAL




Ash




Clay

Iron
METAL



Inox



Black

Om
WOOD



Oak




Wenge


Techno
PMMA



Fluor



Fuxia




Ice

S503



Dimensions



General Specifications

MOTORS
3 engines

DIMENSIONS (LENGTH X WIDTH)
192 x 83 cm

HEIGHT (MIN. AND MAX.)
63/94 cm

WEIGHT
105 kg

THICKNESS
10 cm

DENSITY
30 kg/m³

LIFTING CAPACITY
175 kg

OPERATION
**Hand controller (included)/
Foot controller (optional)**

POWER SUPPLY
230 V - 50/60 Hz

MAX. INTENSITY
2 A

ELECTRIC MOTOR
6.000 N

CUSHION

100% Vinyl-compound

Permablock 3 advanced vinyl protection (against germs, abrasion), UV resistant, superior abrasion 300K Grey colour

Permaguard advanced vinyl protection (against germs and stains), UV resistant, abrasion 100K White colour